
Modulare ^{CVC/PICC} Stabilization Device Instructions for Use

Description: The Modulare Stabilization Device is intended for the securement of catheter hubs with integral suture tabs. Each half of the paired device is comprised of a capture dome on an adhesive pad. The dome and pad assembly is attached to each side of the catheter hub and adhered to the patient’s skin.








Indications for Use: The Modulare device is a stabilization device for suture winged catheters. **Discard after single-patient use.**

Contraindications: Known tape or adhesive allergies.

Definitions

Warning: Could result in death or injury to the patient.
Caution: May result in minor injury to the patient
Note: May result in property damage or advises on the proper use of the product.
Important Technique: Key instruction aiding the successful use of the product.

Warnings, Cautions, and Notes

Symbols	
REF	Catalog Number
LOT	Lot Number
	Do not reuse. For single use only
	Attention: See instructions for Use
	Sterilized by Gamma Radiation. The product is sterile if the product is not open or damaged
QTY	Quantity, Contents
	EU-Conformity Symbol
	Use By
	Manufacturer
	European Representative
	Contains no natural rubber latex
R	Caution: Federal (USA) law restricts this device to sale by or on the order of a Physician

Warning: Read instructions for use prior to using this device.
Warning: Procedure must be performed by trained personnel with knowledge of anatomical landmarks, safe technique and potential complications.
Warning: The Modulare Stabilization Device is provided sterile and is intended for single use only. Do not resterilize. Opened packages and unused product should be discarded.
Warning: Reuse, reprocessing or resterilization may compromise the structural integrity of the devices causing device failure and/or lead to patient infection or cross infection which in turn may result in patient injury, illness or death.
Warning: Before use, inspect all packaging to verify that no damage has occurred during shipping or handling. Carefully inspect the packaging to ensure sterility has not been compromised and verify that the product has not expired.
Caution: Use standard sterile technique when opening product.
Caution: Do not use the Modulare Stabilization Device where loss of adherence could occur, such as with a confused patient, diaphoretic, or non-adherent skin, or when the access device is not monitored daily.
Caution: Observe universal blood and body fluid precautions and infection control procedures, during application and removal of the Modulare device.
Caution: Suture the Modulare Stabilization Device pad to skin if desired or deemed necessary.
Caution: Avoid contact with alcohol or acetone, both can weaken bonding of components and the Modulare Stabilization Device pad adherence.
Caution: Minimize catheter/tube manipulation during application and removal of Modulare Stabilization Device.
Caution: Always apply adhesive strip to catheter near insertion site during device dressing change.
Caution: Always prepare targeted stabilization site with alcohol to degrease skin. Then apply skin prep for enhanced adherence and skin protection. Allow to dry completely (30 seconds).
Caution: Always place sterile transparent dressing over the entire Modulare Stabilization Device and catheter entry site.
Caution: The Modulare Stabilization Device should be replaced at least every 7 days.
Caution: Federal (U.S.A) law restricts this device to sale by or on the order of a physician.

Note: The Modulare Stabilization Device can be adhered after X-ray confirmation of proper catheter/tube position.
Note: After use, this product should be considered a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal regulations.
Technique for use of the Modulare device
Important Technique: Prepare targeted stabilization site with alcohol to degrease skin. Then apply skin prep for enhanced adherence and skin protection. Allow to dry completely (30 seconds).

- Remove the first section of paper backing from one of the Modulare pads exposing the adhesive under the dome.
- Capture one of the catheter’s suture tabs by inserting the Modulare dome post through the suture hole and adhere the exposed adhesive to patient’s skin.
- With the exposed part of the pad adhered, pull back the remaining paper backing while applying downward pressure on the pad to adhere the remainder of the pad.
- Repeat the steps above for the other catheter suture tab. Verify that both dome posts are engaged in the suture tab holes.
- Apply adhesive strip over catheter at or near insertion site.
- Apply sterile transparent dressing over the catheter entry site and the entire Modulare Stabilization Device.
- Carefully remove the Modulare Securement Device with alcohol to loosen adhesive.

Storage: Products should be stored in a clean, cool, dry area away from chemical fumes.

Shelf Life: A “Use By” date is referenced on the label.

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